medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
 Unit 4/13-21 Hallmark Street
 Pendle Hill NSW 2145
 Australia
 T +61 2 9890 8696
 F +61 2 9890 8439
 sales@medi.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: +32-11 24 25 60
 F: +32-11 24 25 64
 info@medi-belgium.be
 www.medi-belgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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Medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
 06510 Çankaya Ankara
 Turkey
 T: +90 312 435 20 26
 F: +90 312 434 17 67
 info@medi-turk.com
 www.medi-turk.com

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 0EL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

4 051526 912366

E013644/04.2025

circaid® juxtafit® essentials upper leg with knee inelastic compression system

Gebrauchsanweisung. Instructions for use. Instructions d'utilisation.
 Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso.
 Bruksveiledning. Brugsanvisning. Návod k použití. دليل الاستخدام. Инструкция
 по использованию. Instrukcja zakładania. Instrukcja użytkowania. Kasutusjuhend.
 Bruksanvisning. Käyttöohje. 使用説明書.

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Important notes

This medical device **[MD]** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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Nylon, Polyurethan, Elastan, Polyethylen

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ENGLISH / ENGLISH

circaid® juxtafit® essentials upper leg with knee

Intended purpose

This compression product serves to compress the leg in patients with venous and lymphatic diseases .

Indications

- Lymphoedema
- Lipoedema
- Chronic venous insufficiency
 - Varicose veins
 - Lipodermatosclerosis
 - Venous ulcer
- Deep vein thrombosis/thromboprophylaxis
- Acute deep vein thrombosis
- Post-thrombotic syndrome
- Post sclerotherapy
- Dependent oedema

Contraindications

- Advanced peripheral arterial disease (if one of these parameters is present: ABPI < 0.5; ankle arterial pressure < 60 mmHg; toe pressure < 30 mmHg or TcPO₂ < 20 mmHg on dorsum of foot). If non-elastic materials are used, the patient can still try out compression clothing if they have ankle arterial pressure between 50 and 60 mmHg, under close clinical supervision.
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Untreated infection of the leg and/or foot
- Any circumstances under which an increase in venous or lymphatic reflux is undesirable.
- Suspected or known untreated acute deep vein thrombosis

Risks/side effects in the event of

- Slight or moderate peripheral arterial disease
- Decreased sensation – the patient must be able to feel the degree of compression
- Sensitivity to one of the used materials
- Use in children and patients dependant on assistance – patients must be able to feel and communicate the applied degree of compression

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Recommendations for wearing

Always make sure that the correct prescribed pressure is applied for compression. At night, slightly loosen the straps of the compression product. Every vertical BPS line on the straps can be at a different spot on the leg according to the size and shape of the leg. They do not have to be aligned to the front centre of the leg or to another defined area of the leg. In case of pain, remove the compression product immediately and consult with your physician.

If applied correctly, the straps overlap slightly so that there are no gaps. Only tighten the straps in the intended places. The hook-and-loop fastener must never come into contact with skin!

Fitting instructions

- Pull up the circaid understocking until it covers the entire leg up to the groin. The foot seam should point outwards. Make sure there are no creases in the understocking (Fig. 1).

Note: If you also use a circaid juxtafit for the lower leg (not included), put this on first. Please follow the instructions for use enclosed with this product.

- Position the compression product just below the groin (Fig. 2). The black side of the compression product must face towards the leg. Make sure that the raised area of the compression product is positioned on the outer side of the leg (Fig. 3).
- Unroll the two top straps and fasten them loosely on the opposite side to fix the compression product in the correct

position (Fig. 4). This makes adjusting the compression product easier.

- Unroll the two lowest straps and fasten the lowest strap on the opposite side of the compression product with a tight, but comfortable, degree of compression (Fig. 5).
- Before fastening the second strap from the bottom, loosen the third strap from the bottom. Now fasten the second strap. Repeat this process with the remaining straps (Figs. 6+7).
- The compression product should fit snugly and smoothly around the leg. The applied degree of compression should be tight but comfortable. Read the section "Using the Built-In Pressure Card" for information on adjusting the degree of compression. Make sure that all straps have been fastened in an alternating order.
- Check the compression product for gaps and creases after all straps have been fastened. If necessary, adjust them (Fig. 8).
- Putting on the knee section: Fasten the knee cover to the compression product in such a way that it encloses your knee tightly and comfortably. Make sure that the knee section perfectly fits your knee (Fig. 9). You have now put on the compression product fully (Fig. 10).

In order to maintain the range of compression throughout the day, simply adjust the individual straps gradually over the course of the day.

Recommendation for the combined use of upper and lower leg compression products with knee: If you combine a circaid juxtafit essentials for the lower leg with a circaid juxtafit essentials for the upper leg with knee, make sure that the upper leg compression product with knee overlaps with the lower leg compression product.

Using the Built-In Pressure card

The set compression value can vary due to movement, specific material characteristics and individual body anatomy.

- Pick up the BPS card. Use the edge of the card that corresponds to the circumference of your ankle (CB). Identify the degree of compression prescribed for you on the colour coded card.
- Start with the lowest strap on the compression product. Place the black triangle on the BPS card against one of the BPS lines on the lowest strap.
- Check where the second BPS line of the strap meets the compression areas on the BPS card (Fig. A).
- If the BPS line is below (Fig. B) or above (Fig. C) the desired compression area, loosen and fasten the strap again until reaching the desired compression (Fig. A). If there is too much elasticity, completely unfasten the strap before tightening again. Make sure that all straps are adjusted to the same level of compression to achieve a decreasing rate of compression from distal to proximal. In order to maintain the prescribed range of compression throughout the day, check the distances between the BPS lines with the BPS card, if necessary. In the event of possible changes, please make adjustments back to the desired compression range.

Please make sure to always cover the entire limb with the medical compression product.

Adjusting the strap length

If a strap is too long, you can shorten it individually. However, a strap must never be cut shorter than is marked by the BPS lines. We recommend that you shorten in small steps to prevent cutting off too much at once. A strap is too long and

should be shortened if the straps can be pulled to the back of the leg or if there is excessive wrinkling on the back of the leg (Fig. 14). In order to shorten a strap, first determine by how much it must be shortened. Mark the point you have determined. Then place the compression product flat on a table or similar surface and cut in a straight line along the strap from the marked point (Fig. 15). Then fasten the edge of the strap with one of the hook and loop fasteners included with the product (Fig. 16). Never cut the compression product while it is being worn.

Care instructions

This compressive product can be washed on a delicate cycle and dried in a tumble dryer set to a low heat. Please roll up the straps for laundry, the hook-and-loop-fasteners must not be open. We recommend using a mesh laundry bag. To extend the life of the product, hand wash and drip dry.

-  Delicate wash cycle
-  Do not bleach
-  Tumble dry at low temperature
-  Do not iron
-  Do not dry-clean
-  Do not wring

Storage

Please ensure that the product is dry before storing it and store it out of direct sunlight.

Product life

6 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition

Nylon, polyurethane, elastane, polyethylene

Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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