sigvaris

Compreflex® Standard Calf

Compreflex® Transition Calf

Compreflex® Complete Calf

Coolflex™ Standard Calf

Gebrauchsanweisung	2	Petunjuk Penggunaan	46
Instructions for Use	4	使用指南	48
Istruzioni per l'uso	6	사용 설명서	50
Mode d'emploi	9	使用方法	51
Brugsanvisning	11	คำแนะนํการใช้งาน	53
Käyttöohje	13	Пайдалану жөнiндегi	
Gebruiksaanwijzing	15	нұсқаулық	55
Bruksanvisning	18	उपयोग के लिए निर्देश	58
Instrucciones de uso	20	Arahan Penggunaan	59
Instruções de utilização	22	62	تعليمات الاستخدام
Инструкция по применению	25	64	הנחיות לשימוש
Instrukcja użytkowania	27		
Návod na použitie	30		
Navodila za uporabo	32		
Οδηγίες χρήσης	34		
Kullanım talimatı	36		
Инструкции за употреба	39		
Instrucțiuni de utilizare	41		
Naudojimo instrukcija	44		

Produktbeschreibung

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English

Product description

The Compreflex® Standard Calf is an inelastic product that provides precise compression ranges of 21–46 mmHg thanks to the use of the Accutabs. The product is combined with a pair of Basic Liners to protect the skin.

The Compreflex® Transition Calf is an inelastic product that provides precise compression ranges of 21–46 mmHg thanks to the use of the Accutabs. The product is combined with a pair of Transition Liners to protect the skin, which apply 10–15 mmHg of compression pressure to the foot and ankle.

The Compreflex® Complete Calf is an inelastic product that provides precise compression ranges of 21–46 mmHg thanks to the use of the Accutabs. The product is combined with a pair of Complete Liners to protect the skin, which apply 10–15 mmHg of compression pressure to the foot, ankle and calf.

The Coolflex™ Standard Calf is an innovative compression device that features a breathable, lightweight material on the one side and a hook-and-eye fastener on the other. Its hook-and-eye strip makes it ideal for users with limited mobility. It can be applied and removed without assistance. In addition, the product features an integrated pressure adjustment system to adjust the compression ranges between 20–50 mmHg. The product is combined with a pair of Basic Liners to protect the skin.

Intended purpose

These inelastic wraps are compression therapy devices used on the calf. Compression therapy improves venous and lymphatic return, which reduces venous and lymphatic edema. It is used in the initial decongestion phase of lymphatic edema and pronounced venous edema as an alternative to wrapping with bandages. This medical device is intended for multiple use and is not sterile.

Effect of medical compression

Compression products apply a carefully controlled pressure on tissue and veins. This improves the blood return in the veins and increases the lymphatic flow. At the same time, it improves microcirculation and therefore the supply to the tissue.

Target group

The intended target group includes patients who have been prescribed a compression product for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions who apply the products to patients.

Indications and contraindications

Indications

- Edema caused by chronic venous disease (CVD) due to obstruction or valvular insufficiency: treatment and/or prevention
- CVD with venous edema, skin changes, open and closed ulcers
- Lymphedema
- Lipedema
- Edema due to immobilization
- Non-specific edema or edema after surgery

Contraindications

- Advanced peripheral arterial occlusive disease PAOD (ABPI <0.6, ankle blood pressure <60 mmHg, toe blood pressure <30 mmHg, or TcPO₂ <20 mmHg instep)
- Severe cardiac insufficiency (NYHA III–IV)
- Septic phlebitis
- Phlegmasia coerulea dolens
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Compression of an existing arterial bypass (epifascial)

Relative contraindications that require a particular balancing of benefits and risks:

- Allergy against compression material
- Inflammatory processes (compression products to reduce inflammation, pain or edema in leg erysipelas or cellulitis should only be used in combination with antibacterial treatment)
- PAOD with ABPI between 0.6 and 0.9 and/or ankle pressure between 60 and 90 mmHg
- Skin diseases with massive exudation
- Primary chronic polyarthritis

Risks and side effects

Inelastic wraps can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Itching, scaling and signs of inflammation may occur under the compression aids in cases of sensitive skin, so adequate skin care is advised. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing, immediately consult your physician. In this case, the compression aid should no longer be worn until otherwise prescribed by a physician.

Please note: Consult your physician. The precautions and contraindications listed here are guidelines only and are no substitute for medical advice.

Special instructions

- Protect the compression product from ointments, oils and greases to preserve its service life.
- Compression products must not be worn over open wounds. Wounds should be covered with appropriate dressings.
- Failure to comply with this instruction may reduce the safety and effectiveness of the product.
- All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.

- Compression is increased when liners with compression pressure are used under a wrap.
- It should be ensured that the pressure of the compression product decreases (degressive pressure gradient) from distal to proximal (away from the body to close to the body).
- Make sure that the liner does not crease.

Application

Applying the Compreflex® Standard Calf/Compreflex® Transition Calf/Compreflex® Complete Calf

1. Pull the liner (included) evenly over your heel and distribute the knitted section evenly over the leg. The optimum fit of the liner is achieved when it finishes 2 finger widths below the back of the knee.
2. Before applying, unfold each hook fastener on the product to its full length, roll the ends back and secure them there so that they cannot get tangled. Slide the foot into the opening between the elastic insert and the product and pull the product up to just below the knee, but no further than the liner. The textile label is located on the top at the back of the leg and faces outward.
3. Carefully position the two lowest hook fasteners over the leg and fasten them without applying compression. The hook fasteners should overlap as little as possible without leaving any gap between them. Press the ends of the fasteners firmly to the sides. Continue positioning the remaining hook fasteners in the same way and fix them in place.
4. Apply the Accutabs: your prescribing physician should tell you which compression range is best for you. Apply the Accutabs accordingly and then close the hook fasteners as directed (refer to the Accutab® instructions for use).

Removing the Compreflex® Standard Calf/Compreflex® Transition Calf/Compreflex® Complete Calf

To remove, loosen each hook fastener, roll their ends backwards and secure them there so that they cannot get tangled. Remove the product from the leg and take off the liner.

Applying the Coolflex™ Standard Calf

1. Pull the liner (included) evenly over your heel and distribute the knitted section evenly over the leg. The optimum fit of the liner is achieved when it finishes 2 finger widths below the back of the knee.
2. Slide the foot into the opening between the elastic insert and the product and pull the product up to just below the knee, but no further than the liner. The textile label is located on the top at the back of the leg and faces outward.
3. Make sure both the top label and the white marker at the bottom are aligned with the shin bone.
4. Starting at the bottom at the ankle, hook the hook-and-eye fasteners into the mesh material and continue upward to below the knee.
5. Tighten all hook-and-eye fasteners to the first (20–30 mmHg pressure) or

second (30–40 mmHg pressure), or third (40–50 mmHg pressure) line as prescribed. The fasteners should be adjusted (pulled tight) while sitting or standing. The selected pressure is the pressure at rest, which may change depending on movement and body position because the product is not very elastic.

6. Check the hook-and-eye fasteners and adjust them as needed, especially in the morning and throughout the day. The applied compression may change the leg circumference and the fasteners should be readjusted as needed to maintain optimal pressure.

Removing the Coolflex™ Standard Calf

To remove, loosen all hook-and-eye fasteners and remove the product from the leg. Take off the liner.

Storage, wearing period and service life

Unless otherwise prescribed by your physician, the compression product can be worn during the day and overnight. However, remove the product once a day to clean the skin and care for any wounds. The product should then be reapplied according to the instructions for use. How often the product needs to be removed for skin checks, hygiene and/or wound care may vary according to the individual needs of patients and the prescriber's recommendations. A liner should be worn under every compression product. Wearing and washing a compression product every day may cause the medical pressure gradient and elasticity of your wrap to diminish. The shelf life date is indicated on the packaging using the symbol of an hourglass. We guarantee maximum effect for the subsequent six months of wearing.

Care and storage

Follow the instructions for care on the packaging and/or the textile label. Close all hook or other fasteners before washing. Wash the product in a mesh laundry bag on a gentle cycle using a mild detergent. Do not use bleach or fabric softener. Dry the product on low/no heat in a dryer or air dry (air drying is recommended to extend wrap life). Store the product at room temperature and in a dry place away from light. Keep away from direct heat. Wash liners before first use.

Disposal

There are no special requirements for disposing of the medical compression product. It can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

List of materials

Compreflex® Standard Calf/Compreflex® Transition Calf/Compreflex® Complete Calf
Polyamide
Polyurethane foam
Cotton
Elastane
Polyester
Made without natural rubber latex

Coolflex™ Standard Calf

Polyester
Polyamide
Cotton
Elastane
Stainless steel
Made without natural rubber latex

Pack contents

Compreflex® Standard Calf

1 Compreflex® Standard Calf
1 pair Basic Liners
1 set Accutab®

Compreflex® Transition Calf

1 Compreflex® Standard Calf
1 pair Transition Liners
1 set Accutab®

Compreflex® Complete Calf

1 Compreflex® Standard Calf
1 pair Complete Liners
1 set Accutab®

Coolflex™ Standard Calf

1 Coolflex™ Standard Calf
1 pair Basic Liners

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2023 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS AG, Gröblistrasse 8, 9014 St. Gallen, Switzerland Phone +41 71 272 40 00
SIGVARIS GmbH, Dr.-Karl-Lenz-Strasse 35, 87700 Memmingen, Germany Phone +49 8331 757 0
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), Italy Phone +39 422 052 750
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, Poland Phone +48 61 43 88 200