by SIGVARIS

Gebrauchs- anweisung	2	Návod k použití	23	Kasutusjuhend	42
Instructions for use	3	Upute za upotrebu	24	Інструкція з користування	44
Istruzioni per l'uso	5	Návod na použitie	26	Petunjuk penggunaan	45
Mode d'emploi	6	Navodila za uporabo	27	使用说明	47
Brugsvejledning	8	Uputstvo za upotrebu	29	사용 설명서	48
Bruksanvisning	10	Οδηγίες χρήσης	31	使用方法	50
Käyttöohjeet	11	Kullanım talimatı	32	คำแนะนำในการใช้งาน	51
Gebruiks- aanwijzing	13	Használati útmutató	31	გამოყენების ინსტრუქცია	52
Bruksanvisning	14	Инструкции за употреба	36	54 تعليمات الاستخدام	
Instrucciones de uso	16			56 הנחיות שימוש	
Instruções de utilização	18	Instrucțiuni de utilizare	37		
Инструкция по применению	19	Lietošanas instrukcija	39		
Instrukcja użytkowania	21	Naudojimo instrukcija	41		

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DELILAH
HIGHLIGHT for men
HIGHLIGHT for women
HIGHLIGHT MARINIÈRE for men
HIGHLIGHT MARINIÈRE for women
MICROFIBER SHADES for men
MICROFIBER SHADES for women
SAMSON
TRAVENO

English

Intended purpose

Support stocking to support venous return. The product is intended for multiple use and for application on the legs (lower leg, thigh, entire leg).

The effects of medical support stockings

The noticeable and controlled compression effect supports the return of blood to the heart, which helps to prevent tired and heavy legs.

Target group

Support stockings are designed for people who have either a standing or sedentary job, who travel regularly (car, train, and plane), and who want to take preventive action to protect their health.

Indications & contraindications

Indications

Prevention of tired and heavy legs during long periods of sitting or standing, even without the presence of varicose veins, prevention of latent swelling, subsequent treatment, and prevention following successful varicose vein treatment.

Contraindications

Medical support stockings must not be worn in the following cases:

- Advanced peripheral arterial occlusive disease (ABPI < 0.6, ankle blood pressure < 60 mmHg, toe blood pressure < 30 mmHg, or TcPO₂ < 20 mmHg instep)
- Severe cardiac insufficiency (NYHA III - IV)
- Septic phlebitis
- Phlegmasia coerulea dolens

- Compression of an existing arterial bypass (epifascial)

Relative contraindications that require a particular balancing of benefits and risks:

- Allergy against compression material
- Inflammatory processes (compression products to reduce inflammation, pain or edema in leg erysipelas or cellulitis should only be used in combination with antibacterial treatment)
- PAOD with ABPI between 0.6 and 0.9 and/or ankle pressure between 60 and 90 mmHg
- Skin diseases with massive exudation
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g., in diabetes mellitus)
- Primary chronic polyarthritis

Risks and side effects

Medical support stockings can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. If problems should arise, the support stockings should no longer be worn.

Itching, scaling, and inflammation may occur under the support stockings in cases of sensitive skin, so adequate skin care under the compression aids is advised.

Special instructions

Do not cut loose ends on the insides of the stocking.

Protect the support stockings from ointments, oils, and grease to preserve their service life.

Support stockings must not be worn over open wounds.

For products with grip tops: regularly ventilate the skin.

Noncompliance with this information may limit safety and effectiveness of the product. All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.

If the product is damaged, its medical effectiveness is no longer guaranteed.

Application

We recommend wearing special gloves to make it easier to put the stockings on and take them off*. Sharp fingernails and rings can cause defects.

* **Note:** gloves should not be worn for MICROFIBER SHADES for women and for men or for HIGHLIGHT MARINIÈRE for women and for men.

Donning

1. The best time is when rising in the morning, when the swelling in your legs has subsided. For easier application, make sure that your legs are dry.
2. Slip the stocking evenly over your heel. Do not pull, but spread the fabric evenly over the leg using the gloves or your hand. While doing so, check for a proper fit at the toe and heel.
3. To protect the seam and the grip top, carefully grasp the knitted fabric and not the seam when donning the support stocking.
4. If it needs repositioning, work the fabric upwards towards the knee with the palm of your hand and spread evenly throughout.

Correct fit

- Knee high stockings:
2 finger widths below the back of the knee
- Thigh stockings:
3 finger widths below the buttocks/crotch

Doffing

1. Take hold of the stockings at the upper edge with gloves or the palm of your hand and pull downwards towards the heel.
2. Work the fabric over the heel with the flat of the hand.

For advice on donning and doffing our support stockings, we recommend watching the videos available at www.sigvaris.com/videos.

Storage, wearing period and service life

We recommend wearing the support stockings daily from morning to night. Wearing and washing your support stocking every day may cause the medical pressure gradient and elasticity of your support stocking to diminish.

We guarantee maximum effect for six months of wearing. The shelf life date is indicated on the packaging using the symbol of an hourglass. After this date, the stockings can be used for the wearing time recommended above.

Care and storage

Follow the instructions for care on the packaging and/or the textile labels. Wash inside a laundry net. The stocking must be stored at room temperature and in a dry place away from light. Keep away from direct heat.

Correct care

Applicable for MICROFIBER SHADES for men / MICROFIBER SHADES for women / TRAVENO

Applicable for HIGHLIGHT MARINIÈRE for men / HIGHLIGHT MARINIÈRE for women / HIGHLIGHT for men / HIGHLIGHT for women / DELILAH / SAMSON

Disposal

There are no special requirements for the disposal of medical support stockings. These can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

Material abbreviations:

CO	cotton
EL	elastane
ED	elastodiene
PA	polyamide
PP	polypropylene
PES	polyester
CMD	modal

Applies to the following products

DELILAH
HIGHLIGHT for men
HIGHLIGHT for women
HIGHLIGHT MARINIÈRE for men
HIGHLIGHT MARINIÈRE for women
MICROFIBER SHADES for men
MICROFIBER SHADES for women
SAMSON
TRAVENO

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Attiecas uz MICROFIBER SHADES for men / MICROFIBER SHADES for women / TRAVENO

Attiecas uz HIGHLIGHT MARINIÈRE for men / HIGHLIGHT MARINIÈRE for women / HIGHLIGHT for men / HIGHLIGHT for women / DELILAH / SAMSON

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DELILAH
HIGHLIGHT for men
HIGHLIGHT for women
HIGHLIGHT MARINIÈRE for men

HIGHLIGHT MARINIÈRE for women
MICROFIBER SHADES for men
MICROFIBER SHADES for women
SAMSON
TRAVENO

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2024 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS GmbH, Dr.-Karl-Lenz-Strasse 35, 87700 Memmingen, **Germany**
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), **Italy**
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, **Poland**
SIGVARIS (Australia) Pty. Ltd., 8/53 Stanley Street, West Melbourne VIC 3003, **Australia**